"Noticeable yet natural"











No need for repetitive and costly fillers

Designed exclusively for the lip

- · Is not an injectable filler
- · Contoured to match the shape of your lip
- · Soft, solid silicone-based device
- · Is permanent yet can easily be removed
- Available in various sizes

The PermaLip™ lip implant is available in a variety of sizes to accommodate individual needs ... ranging from slight accent to full, voluptuous lips.



Clinic:

Surgisol.com





"The safe, definitive solution for permanent lip enhancement"

About Lip Enhancement

One of the greatest challenges in plastic and cosmetic surgery of the perioral region has been the search for an exceptional solution to lip enhancement.

Past methodologies, both temporary and permanent, have met with disappointing results for both physician and patient.

Temporary solutions typically involve the use of commercially produced fillers or autologuous fat, which can be unpredictable with respect to how long they last, and require treatments that are painful, repetitive, and costly.

Until now; permanent solutions relied on implants that may, due to their surface characteristics, promote tissue ingrowth. This may result in a restriction of movement and an unnatural appearance.

PermaLip[™] has now addressed these issues to give you that noticeable, yet natural look. It is also the only implant of its type that is both CE and FDA approved.

"The alternative to injectable fillers"

The Solution

PermaLip[™] is made of a soft, solid silicone elastomer using an injection moulding process to produce a naturally contoured implant with a smooth, non-porous surface.

The implant cannot rupture or deflate, and may be removed at any time.

PermaLip™ is the answer.

"Lip enhancement the way it was always meant to be"



The Procedure

PermaLip™ lip enhancement is a simple and quick procedure to perform. The mouth and lips are locally anaesthetised; then two small incisions are made in the corners of the mouth.

The PermaLip™ implant is then passed through a tunnel in the lip and centered. The process is repeated for the other lip and the incisions are closed using an absorbable suture.

The entire process is designed to take less than 15 minutes and normal activity may be resumed the following day.

Reversible

Although PermaLip™ is permanent, it is also designed to be easily reversible should a person wish to return to their original lip appearance. Removal of the implant(s) typically takes less than 10 minutes.

"Enjoy fuller, more youthful lips for years to come"